

# Nexxo-Prep Plasmide mini

DNA extraction, by spin-column system, for the isolation, in less than 15 minutes, of up to 20 µg plasmid DNA from 0.5 – 2.0 ml bacterial culture.

Eluted plasmid DNA is ready for down-stream applications and can directly be used for PCR, sequencing, enzymatic digestion, cloning, etc., or can be stored for future use.

## I. Kit components

	10 preps	250 preps
<b>Tampon d'éluion P</b> (Elution Buffer P)	2 ml	30 ml
<b>Solution de lavage SLB</b> (Wash Solution SLB)	10 ml (ready-to-use)	40 ml (final volume: 200 ml)
<b>Tampon de resuspension</b> (Resuspension Buffer)	2 x 2 ml	70 ml
<b>Tampon de lyse L</b> (Lysis Buffer L)	2 x 2 ml	70 ml
<b>Tampon de neutralisation N</b> (Neutralization Buffer N)	2 x 2 ml	70 ml
<b>Filtres de centrifugation</b> (Spin Filter)	10	5 x 50
<b>Tubes receveurs 1,5 ml</b> (1,5 ml Receiver Tubes)	10	5 x 50
<b>Tubes receveurs 2,0 ml</b> (2,0 ml Receiver Tubes)	10	5 x 50
<b>User guide</b>	1	1
<b>Art. No.</b>	2036.10	2036.250

### Required material and equipment not included in this kit

- Ethanol >96 %
- Reaction tubes (1.5 ml / 2.0 ml)
- Microcentrifuge (12000 - 16000 x g)
- Pipettes with corresponding tips
- Disposable gloves
- Lysozyme (only for protocol 3)

Some components are delivered in concentrated form and have to be diluted appropriately (see chapter « Reagents and buffer solutions preparation », page 2).

## II. Storage and stability

---

All kit components should be stored at room temperature (15-30 °C).

Ethanol is a volatile compound. Keep **Wash Solution SLB** tightly closed.

Check solutions for absence of precipitates before use. Redissolve precipitates by heating (< 37 °C).

*Note: all kit components are stable for at least 12 months.*

## III. Reagents and buffer solutions preparation

---

### 1. Kit 10 extractions:

All components are delivered ready-to-use.

### 2. Kit 250 extractions:

Add 160 ml of >96 % ethanol to the **Wash Solution SLB** and store the bottle tightly closed.

## IV. Protocol 1: plasmid extraction from 0.5 - 2.0 ml bacterial culture (ex.: *E. coli*)

### **Before starting**

- Ensure that reagents/buffers preparation has been done (see chapter "Reagents and buffer solutions preparation", page 2).
- Ensure that all components are at room temperature and check solutions for absence of precipitates before use. If necessary resuspend precipitates by heating (< 37 °C).

*Note: To prevent contamination, use new pipet tip for each pipetting step.*

<b>1</b>	<b>Preparation of bacterial culture</b>
	<ul style="list-style-type: none"> <li>• Inoculate 1 – 5 ml of LB broth (supplemented with the appropriate selective antibiotic), with a single colony</li> <li>• Incubate 12-16 hours at 37 °C with rapid shaking</li> </ul> <p><i>Note: incubation period should not be longer than 16 hours.</i></p>

<b>2</b>	<b>Bacterial cells harvesting</b>
	<ul style="list-style-type: none"> <li>• Transfer 0.5 - 2.0 ml of bacterial culture in a 1.5 or 2.0 ml reaction tube</li> <li>• Centrifuge 1 min. at max. speed (12000 - 16000 x g)</li> <li>• Discard the whole supernatant, and keep only the pellet</li> </ul>

<b>3</b>	<b>Resuspension of bacterial cells</b>
	<ul style="list-style-type: none"> <li>• Add 250 µl of <b>Resuspension Buffer</b></li> <li>• Vortex until pellet is completely resuspended (any pellets or cell clumps should remain)</li> </ul>

<b>4</b>	<b>Cell lysis and neutralization</b>
	<ul style="list-style-type: none"> <li>• Add 250 µl of <b>Lysis Buffer L</b></li> <li>• Close the tube and mix <u>gently</u> by inversion (5 inversions, do not vortex)</li> </ul>
	<hr/> <p>Caution: lysis step (↑) should not take more than 5 minutes.</p> <hr/> <ul style="list-style-type: none"> <li>• Add 250 µl of <b>Neutralization Buffer N</b></li> <li>• Mix by inversion (mix <u>gently</u> but effectively, 4 – 6 inversions)</li> <li>• Centrifuge 5 min. at max. speed (12000 - 16000 x g)</li> </ul>

Steps 5 to 7 →

<b>5</b>	<b>Plasmid DNA adsorption to Spin Filter</b>
	<ul style="list-style-type: none"> <li>• Insert a <b>Spin Filter</b> into a new <b>2,0 ml Receiver Tube</b></li> <li>• Transfer the clarified supernatant into the <b>Spin Filter</b></li> <li>• Incubate 1 min. at room temperature</li> <li>• Centrifuge 1 min. at 11000 x g</li> <li>• Discard the flow-through and put the <b>Spin filter</b> back into the <b>2.0 ml Receiver tube</b></li> </ul>

<b>6</b>	<b>Plasmid DNA washing</b>
	<ul style="list-style-type: none"> <li>• Add 750 µl of <b>Wash Solution SLB</b> to the <b>Spin Filter</b></li> <li>• Centrifuge 1 min. at 11000 x g</li> <li>• Discard the flow-through and put the <b>Spin filter</b> back into the <b>2.0 ml Receiver tube</b></li> <li>• Centrifuge 3 min. at max. speed (12000 - 16000 x g), to remove remaining ethanol from the <b>Spin Filter</b></li> </ul>

<b>7</b>	<b>Plasmid DNA elution</b>
	<ul style="list-style-type: none"> <li>• Insert the <b>Spin filter</b> into a new <b>1.5 ml Receiver tube</b></li> <li>• Add 50-100 µl of <b>Elution Buffer P</b> onto the center of the <b>Spin Filter</b></li> <li>• Incubate 1 min. at room temperature</li> <li>• Centrifuge 1 min. at 11000 x g, to elute plasmid DNA</li> </ul> <p><i>Note: depending on desired yield and concentration, DNA can be eluted with more or less elution buffer.</i></p> <p><i>Higher yield is reached when eluting twice (2 x 100 µl for instance).</i></p>

Eluted plasmid DNA is ready for down-stream applications and can directly be used for PCR, sequencing, enzymatic digestion, cloning, etc., or can be stored at 2 °C to 8 °C.

Storing pDNA at -15 to -25°C might be damaging (shearing of pDNA). Repeated freeze-thaw cycles enhance the risk of pDNA degradation.

*Note: **Elution buffer P** contains EDTA. Elution can also be achieved with ddH<sub>2</sub>O. Without buffer pDNA may deteriorate. If pDNA is eluted with ddH<sub>2</sub>O it is recommended to store pDNA at -20 °C.*

## V. Protocol 2: extraction of low copy plasmids and cosmids, from up to 10 ml bacterial culture

*Note: this protocol needs a higher volume of buffers. Using this protocol reduces the number of extractions.*

### **Before starting**

- Ensure that reagents/buffers preparation has been done (see chapter “Reagents and buffer solutions preparation”, page 2).
- Ensure that all components are at room temperature and check solutions for absence of precipitates before use. If necessary resuspend precipitates by heating (< 37 °C).
- For plasmids or cosmids with more than 10 kb, the last step needs **Elution Buffer P** (or ddH<sub>2</sub>O, as the case may be) heated at 70 °C. Warm up in due time required volume of **Elution Buffer P** in a 2 ml tube (tube not supplied).

*Note: To prevent contamination, use new pipet tip for each pipetting step.*

1	Preparation of bacterial culture
	<ul style="list-style-type: none"><li>• Inoculate 1 – 10 ml of LB broth (supplemented with the appropriate selective antibiotic), with a single colony</li><li>• Incubate 12-16 hours at 37 °C with rapid shaking</li></ul> <p><i>Note: incubation period should not be longer than 16 hours.</i></p>

3	Resuspension of bacterial cells
	<ul style="list-style-type: none"><li>• Add 500 µl of <b>Resuspension Buffer</b></li><li>• Vortex until pellet is completely resuspended (any pellets or cell clusters should remain)</li><li>• Transfer the whole solution in a 2.0 ml reaction tube (not supplied)</li></ul>

2	Bacterial cells harvesting
	<ul style="list-style-type: none"><li>• Transfer 1 - 10 ml of bacterial culture in a 15 ml Falcon reaction tube</li><li>• Centrifuge 10 min. at 5400 x g and 4 °C</li><li>• Discard the whole supernatant, and keep only the pellet</li></ul>

Steps 4 to 7 →

<b>4</b>	<b>Cell lysis and neutralization</b>
	<ul style="list-style-type: none"> <li>• Add 500 µl of <b>Lysis Buffer L</b></li> <li>• Close the tube and mix <u>gently</u> by inversion (5 inversions, do not vortex)</li> </ul>
	<p>Caution: lysis step (↑) should not take more than 5 minutes.</p> <ul style="list-style-type: none"> <li>• Add 500 µl of <b>Neutralization Buffer N</b></li> <li>• Mix by inversion (mix <u>gently</u> but effectively, 4 – 6 inversions)</li> <li>• Centrifuge 5 min. at max. speed (12000 - 16000 x g)</li> </ul>

<b>5</b>	<b>Plasmid DNA adsorption to Spin Filter</b>
	<ul style="list-style-type: none"> <li>• Insert a <b>Spin Filter</b> into a new <b>2,0 ml Receiver Tube</b></li> <li>• Transfer 750 µl of clarified supernatant, from the previous step, into the <b>Spin Filter</b></li> <li>• Incubate 1 min. at room temperature</li> <li>• Centrifuge 1 min. at 11000 x g</li> <li>• Discard the flow-through and put the <b>Spin filter</b> back into the <b>2.0 ml Receiver tube</b></li> <li>• Transfer the leftover clarified supernatant, from step 4, into the <b>Spin Filter</b></li> <li>• Incubate 1 min. at room temperature</li> <li>• Centrifuge 1 min. at 11000 x g</li> <li>• Discard the flow-through and put the <b>Spin filter</b> back into the <b>2.0 ml Receiver tube</b></li> </ul>

<b>6</b>	<b>Plasmid DNA washing</b>
	<ul style="list-style-type: none"> <li>• Add 750 µl of <b>Wash Solution SLB</b> to the <b>Spin Filter</b></li> <li>• Centrifuge 1 min. at 11000 x g</li> <li>• Discard the flow-through and put the <b>Spin filter</b> back into the <b>2.0 ml Receiver tube</b></li> <li>• Centrifuge 3 min. at max. speed (12000 - 16000 x g), to remove remaining ethanol from the <b>Spin Filter</b></li> </ul>

<b>7</b>	<b>Plasmid DNA elution</b>
	<ul style="list-style-type: none"> <li>• Insert the <b>Spin filter</b> into a new <b>1.5 ml Receiver tube</b></li> <li>• Add 50-100 µl of <b>Elution Buffer P (*)</b> onto the center of the <b>Spin Filter</b></li> </ul>
	<p>(*) For plasmids or cosmids with more than 10 kb, use <b>Elution Buffer P</b> heated at 70 °C</p> <ul style="list-style-type: none"> <li>• Incubate 1 min. at room temperature</li> <li>• Centrifuge 1 min. at 11000 x g, to elute plasmid DNA</li> </ul> <p><i>Note: depending on desired yield and concentration, DNA can be eluted with more or less elution buffer.</i></p> <p><i>Higher yield is reached when eluting twice (2 x 100 µl for instance).</i></p>

Eluted plasmid DNA is ready for down-stream applications and can directly be used for PCR, sequencing, enzymatic digestion, cloning, etc., or can be stored at 2 °C to 8 °C.

Storing pDNA at -15 to -25°C might be damaging (shearing of pDNA). Repeated freeze-thaw cycles enhance the risk of pDNA degradation.

*Note: **Elution buffer P** contains EDTA. Elution can also be achieved with ddH<sub>2</sub>O. Without buffer pDNA may deteriorate. If pDNA is eluted with ddH<sub>2</sub>O it is recommended to store pDNA at -20 °C.*

## VI. Protocol 3: plasmid extraction from 0.5 - 2.0 ml Gram positive bacterial culture

### **Before starting**

- Ensure that reagents/buffers preparation has been done (see chapter “Reagents and buffer solutions preparation”, page 2).
- Ensure that all components are at room temperature and check solutions for absence of precipitates before use. If necessary resuspend precipitates by heating (< 37 °C).
- Prepare 10 µl of Lysozyme (10 mg/ml or according to manufacturer’s instructions) for step 3 (Resuspension of bacterial cells)

*Note: To prevent contamination, use new pipet tip for each pipetting step.*

<b>1</b>	<b>Preparation of bacterial culture</b>
	<ul style="list-style-type: none"> <li>• Inoculate 1 – 5 ml of LB broth (supplemented with the appropriate selective antibiotic), with a single colony</li> <li>• Incubate 12-16 hours at 37 °C with rapid shaking</li> </ul> <p><i>Note: incubation period should not be longer than 16 hours.</i></p>

<b>2</b>	<b>Bacterial cells harvesting</b>
	<ul style="list-style-type: none"> <li>• Transfer 0.5 - 2.0 ml of bacterial culture in a 1.5 or 2.0 ml reaction tube</li> <li>• Centrifuge 1 min. at max. speed (12000 - 16000 x g)</li> <li>• Discard the whole supernatant, and keep only the pellet</li> </ul>

<b>3</b>	<b>Resuspension of bacterial cells</b>
	<ul style="list-style-type: none"> <li>• Add 250 µl of <b>Resuspension Buffer</b></li> <li>• Vortex until pellet is completely resuspended (any pellets or cell clusters should remain)</li> <li>• Add 10 µl of Lysozyme (10 mg/ml or according to manufacturer’s instructions)</li> <li>• Vortex</li> <li>• Incubate 10 min. at 37 °C</li> </ul>

Steps 4 to 7 →

<b>4</b>	<b>Cell lysis and neutralization</b>
	<ul style="list-style-type: none"> <li>• Add 250 µl of <b>Lysis Buffer L</b></li> <li>• Close the tube and mix <u>gently</u> by inversion (5 inversions, do not vortex)</li> </ul>
	<p>Caution: lysis step (↑) should not take more than 5 minutes.</p> <ul style="list-style-type: none"> <li>• Add 250 µl of <b>Neutralization Buffer N</b></li> <li>• Mix by inversion (mix <u>gently</u> but effectively, 4 – 6 inversions)</li> <li>• Centrifuge 5 min. at max. speed (12000 - 16000 x g)</li> </ul>

<b>5</b>	<b>Plasmid DNA adsorption to Spin Filter</b>
	<ul style="list-style-type: none"> <li>• Insert a <b>Spin Filter</b> into a new <b>2,0 ml Receiver Tube</b></li> <li>• Transfer the clarified supernatant into the <b>Spin Filter</b></li> <li>• Incubate 1 min. at room temperature</li> <li>• Centrifuge 1 min. at 11000 x g</li> <li>• Discard the flow-through and put the <b>Spin filter</b> back into the <b>2.0 ml Receiver tube</b></li> </ul>

<b>6</b>	<b>Plasmid DNA washing</b>
	<ul style="list-style-type: none"> <li>• Add 750 µl of <b>Wash Solution SLB</b> to the <b>Spin Filter</b></li> <li>• Centrifuge 1 min. at 11000 x g</li> <li>• Discard the flow-through and put the <b>Spin filter</b> back into the <b>2.0 ml Receiver tube</b></li> <li>• Centrifuge 3 min. at max. speed (12000 - 16000 x g), to remove remaining ethanol from the <b>Spin Filter</b></li> </ul>

<b>7</b>	<b>Plasmid DNA elution</b>
	<ul style="list-style-type: none"> <li>• Insert the <b>Spin filter</b> into a new <b>1.5 ml Receiver tube</b></li> <li>• Add 50-100 µl of <b>Elution Buffer P</b> onto the center of the <b>Spin Filter</b></li> <li>• Incubate 1 min. at room temperature</li> <li>• Centrifuge 1 min. at 11000 x g, to elute plasmid DNA</li> </ul> <p><i>Note: depending on desired yield and concentration, DNA can be eluted with more or less elution buffer.</i></p> <p><i>Higher yield is reached when eluting twice (2 x 100 µl for instance).</i></p>

Eluted plasmid DNA is ready for down-stream applications and can directly be used for PCR, sequencing, enzymatic digestion, cloning, etc., or can be stored at 2 °C to 8 °C.

Storing pDNA at -15 to -25°C might be damaging (shearing of pDNA). Repeated freeze-thaw cycles enhance the risk of pDNA degradation.

*Note: **Elution buffer P** contains EDTA. Elution can also be achieved with ddH<sub>2</sub>O. Without buffer pDNA may deteriorate. If pDNA is eluted with ddH<sub>2</sub>O it is recommended to store pDNA at -20 °C.*